



## Thrombocytopenia

# Differential diagnosis of thrombocytopenia

## Thrombocytopenia and automated platelet measurement

Thrombocytopenia is a condition characterised by an abnormally low platelet count – lower than the normal platelet count in adults that ranges from  $150 \times 10^9/L$  to  $450 \times 10^9/L$ . Severe thrombocytopenia, with platelet counts lower than  $20 \times 10^9/L$ , is associated with spontaneous bleeding (not caused by injury). Overlooking a severe thrombocytopenia can have serious consequences for the patient, so obtaining reliable platelet counts is essential for making clinically important decisions with confidence.

Yet obtaining accurate platelet counts, particularly from samples with thrombocytopenic conditions, is not always easy and a challenging task for the laboratory. Sysmex offers a readily available platelet counting solution with PLT-F. You will find more information on this topic in the section 'Challenges in determining an accurate and precise platelet count'.

## Aetiology of thrombocytopenia

Although thrombocytopenia is defined by low platelet concentrations, PLT counts alone do not reveal the underlying causes, which can be inherited or acquired. The causes can be divided into two main categories: decreased bone

marrow production and increased destruction/consumption of platelets in peripheral blood. Often, the clinical question is whether thrombocytopenia is due to bone marrow failure as observed in conditions such as aplastic anaemia (AA) or myelodysplastic syndromes (MDS), or due to increased peripheral destruction/consumption such as in immune thrombocytopenia (ITP), thrombotic thrombocytopenic purpura (TTP) or disseminated intravascular coagulation (DIC). Invasive bone marrow biopsies are usually recommended to investigate the underlying aetiology.

## Differential diagnosis of thrombocytopenia

Differential diagnosis of thrombocytopenia is complex and requires an investigation of the patient's medical history, an evaluation of clinical symptoms, functional platelet tests and an assessment of blood-derived platelet parameters. Historically, some clinicians have used the mean platelet volume (MPV) as a surrogate marker for platelet production because immature platelets tend to be larger than mature platelets. However, the presence of schistocytes, microcytes or other particles with a volume similar to platelets can make the MPV unreliable. Additionally, MPV values are imprecise or impossible to determine, particularly in samples with very low PLT counts, for which information about platelet production is needed most.

The immature platelet fraction (IPF) is a better marker for platelet production and indicates the percentage of immature platelets in relation to the total PLT count. It was described in 1992 by Ault *et al.*, who coined the term 'reticulated platelets' to describe newly released platelets with an elevated RNA content, whose numbers correlated with megakaryocytic activity [1]. IPF is a reproducible parameter that correlates well with the reticulated platelet count obtained from CD61 flow cytometry [2]. Although there is only a partial correlation with the MPV, immature platelets tend to be larger than mature platelets: one study found that 61% of reticulated platelets were in the tertile containing the largest platelets, while 32% and 7% were in the middle and small tertiles, respectively [3]. In addition, immature platelets are more reactive than mature platelets. They contain higher amounts of RNA and are able to produce various proteins typical for active platelets (e. g. GPIIb/IIIa, P-selectin) [3].

### IPF reference ranges

Several studies have established reference ranges for IPF on the Sysmex XE- and XN-Series [4–12]. Generally, there is good consistency among these studies, but in a recent comprehensive study, L van Pelt J *et al.* analysed 12,782 blood samples from Dutch healthy individuals and established reference ranges for IPF in the range of 1.2–8.9% [13]. One study established an IPF reference range for neonates: 0.7–7.9% [14]. However, reference ranges should be always examined for suitability in a given patient population according to the method recommended by the International Federation of Clinical Chemistry and Laboratory Medicine [15].


Thrombocytopenia with an increased IPF may indicate increased destruction in peripheral blood, loss of platelets or a hereditary macrothrombocytopenia [2, 5, 10–12, 16–20, 22].

Thrombocytopenia with a normal or decreased IPF may indicate decreased platelet production in bone marrow [2, 5, 10–12, 16–20].


Several publications reported that IPF obtained from the Sysmex XE- and XN-Series analysers is higher in patients with thrombocytopenia caused by excessive platelet destruction/consumption than in patients with thrombocytopenia caused by decreased platelet production in bone marrow [2, 5, 10–12, 16–20].

### For example:

- Briggs *et al.* followed patients with ITP and TTP, which are both caused by excessive platelet consumption: IPF values were markedly elevated in these two patient groups, while patients under chemotherapy (causing decreased bone marrow production of platelets), and ITP and TTP patients in remission had normal IPF values (Fig. 1) [10].
- Kickler *et al.* reported high IPF values in thrombocytopenic patients with increased platelet destruction, while normal to slightly increased values were observed in patients with decreased platelet production (Table 1) [11].
- Abe *et al.* used a cut-off of 7.7%, resulting in a sensitivity of 86.8% and a specificity of 92.6% for the differential diagnosis of ITP and AA. In addition, IPF was found to be more useful than the mean platelet volume [16].
- Jung *et al.* found that IPF is higher in ITP patients than in AA patients and that a cut-off of 7.3% could be used to distinguish ITP from AA, resulting in a sensitivity of 54.0% and a specificity of 92.2% [12]. (The lower sensitivity compared to the study of Abe *et al.* could be due to different patient cohorts: Patients with acute ITP typically have a high IPF, while patients in remission can also have a normal IPF.)
- Strauss *et al.* studied children with thrombocytopenia and observed that IPF was low in children with platelet production defects, while it was markedly increased in ITP patients, indicating accelerated platelet turnover [18].
- Sakuragi *et al.* found that IPF from the XN-Series had a higher precision and was affected by fewer interferences than IPF from the XE-series. Using a cut-off of 5.8% resulted in a sensitivity of 85.1% and a specificity of 89.3% to distinguish ITP from aplastic thrombocytopenia [20].



**Fig. 1** IPF values in different patient groups. ITP: immune thrombocytopenia; ITP < 50: ITP patients with a PLT count below  $50 \times 10^9/L$ ; TTP: thrombotic thrombocytopenic purpura; Chemo: patients under chemotherapy. Adapted from Briggs *et al.* [10].



**Fig. 2** IPF values in different patient groups. MYH9D: MYH9 disorders; MTPs: macrothrombocytopenias; ITP: immune thrombocytopenia. Adapted from Miyazaki *et al.* [22].

**Table 1** IPF values in different patient groups.

ITP: immune thrombocytopenia; DIC: disseminated intravascular coagulation; AA: aplastic anaemia; PNH: paroxysmal nocturnal haemoglobinuria; NS: not statistically significant compared to control. Adapted from Kickler *et al.* [11].

Subject	Sample size	Mean	p*
<b>Healthy</b>	80	3.1	–
<b>Destruction</b>			
ITP	37	15.0	< .0001
DIC	25	9.5	< .0001
All destruction	62	12.8	< .0001
<b>Suppression</b>			
AA/PNH	3	6.1	.019
Cancer	16	3.8	NS
All suppression	19	4.1	.05

\* Statistical p-value compared to healthy controls

In summary: the IPF parameter provides an assessment of bone marrow platelet production and aids differentiation between thrombocytopenia caused by decreased bone marrow production and thrombocytopenia caused by increased destruction/consumption. It delivers additional information, which might help in differential diagnosis of thrombocytopenia.

### The value of IPF in the differential diagnosis of congenital thrombocytopenia

The IPF can also contribute to the differential diagnosis of suspected congenital thrombocytopenia. Congenital thrombocytopenia is usually suspected in the case of neonatal thrombocytopenia, the onset of bleeding symptoms in childhood, a family history of thrombocytopenia, or when the PLT count is unresponsive to ITP treatment. The MPV is often used for the differential diagnosis of hereditary thrombocytopenia [21] but, as mentioned before, the MPV is affected by interferences and its value is often imprecise or impossible to determine in samples with very low PLT counts.

Several publications describe how IPF can contribute to the differential diagnosis of congenital thrombocytopenia. For example, Miyazaki *et al.* revealed that the IPF was about five times higher in May-Hegglin MYH9 disorders ( $48.6\% \pm 1.9$ ) and approximately double in other macrothrombocytopenia conditions ( $18.4\% \pm 2.1$ ) compared to ITP patients with similar PLT counts ( $9.2\% \pm 0.3$ ) (Fig. 2) [22]. In contrast, patients with Wiskott-Aldrich congenital microthrombocytopenia (WAS) had a lower IPF than would be expected for their level of thrombocytopenia, and the IPF in these patients was lower than in ITP patients [23]. Similar findings were reported in seven children with WAS, who had a lower absolute IPF count than age-matched chronic ITP patients [24].

## Challenges in determining an accurate and precise platelet count

Automated haematology analysers generally deliver an accurate and precise measurement of platelet counts based on the impedance method (PLT-I). However, interfering particles can result in falsely high counts while precision may be limited when patients suffer from severe thrombocytopenia ( $PLT \leq 20 \times 10^9/L$ ) as the low platelet count limits the number of analysed cells. To resolve this, Sysmex XR-Series and XN-Series analysers can perform reflex measurements with alternative flow cytometry methods (PLT-O or PLT-F) in case interfering particles or severe thrombocytopenia are suspected (Table 2).

In XR-Series and XN-Series analysers equipped with PLT-I and PLT-O, a switching algorithm automatically determines whether a PLT-O reflex measurement is required for reporting an accurate PLT count; in analysers equipped with PLT-I and PLT-F, a PLT-F reflex measurement may be required and is performed likewise. Especially if a severe thrombocytopenia is suspected, a more precise measurement is needed to obtain reliable results for confidently making clinically important decisions. Here, PLT-F would be the method of choice.

## Conclusion

PLT counts alone do not reveal the underlying aetiology of thrombocytopenia. The causes of thrombocytopenia can be due to decreased platelet production in bone marrow or an increased destruction/consumption of platelets in peripheral blood. The IPF is a diagnostic parameter that can support treating physicians when determining the cause of thrombocytopenia based on the aetiology of various congenital and acquired thrombocytopenic states, as described in this white paper.

A high IPF suggests consumptive thrombocytopenic disorders or congenital macrothrombocytopenia, and may also indicate an appropriate bone marrow response to thrombocytopenia. In contrast, a low or normal IPF is seen in aplastic states (Table 3). IPF measurements can be ordered for certain patient populations or measured as a reflex test for unknown patients with unclear thrombocytopenia.

**Table 2** Comparison of different platelet counting methods available on the Sysmex XR-Series and XN-Series analysers.

	Impedance count (PLT-I)	Optical count (PLT-O)	Fluorescence count (PLT-F)
<b>Workflow</b>	Default and routine automated method	Often a reflex method	Often a reflex method
<b>Analysis</b>	Part of complete blood count	Part of reticulocyte count	Dedicated PLT count
<b>Precision, accuracy and interferences</b>	<ul style="list-style-type: none"> <li>■ Low precision if <math>PLT &lt; 20 \times 10^9/L</math></li> <li>■ Low accuracy for samples with interferences: containing particles with a volume similar to platelets (reagent crystals, air bubbles, microcytes, RBC fragments)</li> </ul>	<ul style="list-style-type: none"> <li>■ Low precision if <math>PLT &lt; 20 \times 10^9/L</math></li> <li>■ High accuracy when RBC abnormalities are present</li> <li>■ Low accuracy for samples with WBC fragments (apoptosis/ necrosis)</li> </ul>	<ul style="list-style-type: none"> <li>■ High precision down to <math>PLT = 3 \times 10^9/L</math> due to five-fold counting volume</li> <li>■ Virtually no interferences</li> <li>■ Comparable with reference method (CD41/CD61) [25, 26]</li> </ul>
<b>Diagnostic parameters beyond the PLT count</b>	PDW, MPV, PCT, P-LCR	None	IPF, IPF#

**Table 3** Aetiology of thrombocytopenia and associated IPF values. The ranges in the table are based on the literature [13, 22] and provided for guidance only; interpretation of IPF should always occur within the complete clinical context, including clinical symptoms and other laboratory tests.

Acquired		Hereditary
Ineffective platelet production	Increased platelet destruction/consumption	Congenital macrothrombocytopenia
IPF 1.2 – 8.9%	IPF > 8.9%	IPF > 12%
<b>Bone marrow damage</b> <ul style="list-style-type: none"> <li>■ Neoplastic bone marrow infiltration</li> <li>■ Aplastic anaemia caused by chemicals, drugs or infections</li> <li>■ Chronic ITP with apoptotic megakaryocytes</li> </ul>	<b>Immune causes</b> <ul style="list-style-type: none"> <li>■ Immune thrombocytopenia (ITP)</li> <li>■ Heparin-induced thrombocytopenia (HIT) type II</li> </ul>	<b>IPF &gt; 12%</b> <ul style="list-style-type: none"> <li>■ Bernard-Soulier syndrome</li> <li>■ ACTN1-related thrombocytopenia</li> <li>■ <math>\alpha\delta</math>-storage pool disease</li> <li>■ Variant form of Glanzmann thrombasthenia</li> </ul>
<b>Ineffective production</b> <ul style="list-style-type: none"> <li>■ Megaloblastic anaemia</li> </ul>	<b>Non-immune causes</b> <ul style="list-style-type: none"> <li>■ Thrombotic thrombocytopenic purpura (TTP)</li> <li>■ Haemolytic uraemic syndrome (HUS)</li> <li>■ Disseminated intravascular coagulation (DIC)</li> <li>■ HIT type I</li> <li>■ Bleeding</li> </ul>	<b>IPF &gt; 40%</b> <ul style="list-style-type: none"> <li>■ May-Hegglin MYH9 disorders</li> </ul>

## References

- [1] **Ault A et al. (1992):** The significance of platelets with increased RNA content (reticulated platelets). A measure of the rate of thrombopoiesis. *Am J Clin Pathol.* 98(6): 637–46.
- [2] **Pons I et al. (2010):** Correlation between immature platelet fraction and reticulated platelets. Usefulness in the etiology diagnosis of thrombocytopenia. *Eur J Haematol.* 85(2): 158–63.
- [3] **Guthikonda S et al. (2008):** Role of reticulated platelets and platelet size heterogeneity on platelet activity after dual antiplatelet therapy with aspirin and clopidogrel in patients with stable coronary artery disease. *J Am Coll Cardiol.* 52(9): 743–9.
- [4] **Zucker MJ et al. (2006):** Immature platelet fraction as a predictor of platelet recovery following hematopoietic progenitor cell transplantation. *Lab Hematol.* 12(3): 125–30.
- [5] **Cho YG et al. (2007):** Clinical usefulness of the simple technique to diagnose thrombocytopenia using immature platelet fraction. *Korean J Lab Med.* 27(1): 1–6.
- [6] **Pekelharing JM et al. (2010):** Haematology reference intervals for established and novel parameters in healthy adults. *Sysmex J Int.* 20(1): 1–11.
- [7] **Ko YJ et al. (2013):** Establishment of reference interval for immature platelet fraction. *Int J Lab Hematol.* 35(5): 528–33.
- [8] **Sysmex Corporation (2014):** Reference ranges analysis document for XN-Series.
- [9] **Seo A et al. (2015):** Reference intervals for immature platelet fraction and immature platelet count. *Int J Lab Hematol.* 37(1): e1–2.
- [10] **Briggs C et al. (2004):** Assessment of an immature platelet fraction (IPF) in peripheral thrombocytopenia. *Br J Haematol.* 126(1): 93–9.
- [11] **Kickler T et al. (2006):** A clinical evaluation of high fluorescent platelet fraction percentage in thrombocytopenia. *Am J Clin Pathol.* 125(2): 282–7.
- [12] **Jung H et al. (2010):** Immature platelet fraction: establishment of a reference interval and diagnostic measure for thrombocytopenia. *Korean J Lab Med.* 30(5): 451–9.
- [13] **L van Pelt J et al. (2022):** Reference intervals for Sysmex XN hematological parameters as assessed in the Dutch Lifelines cohort. *Clin Chem Lab Med.* 60(6): 907–920.
- [14] **Cremer M et al. (2004):** Immature platelet fraction as novel laboratory parameter predicting the course of neonatal thrombocytopenia. *Br J Haematol.* 144(4): 619–21.
- [15] **Solberg HE. (2004):** The IFCC recommendation on the estimation of reference intervals. The RefVal program. *Clin Chem Lab Med.* 42(7): 710–4.
- [16] **Abe Y et al. (2006):** A simple technique to determine thrombopoiesis level using immature platelet fraction (IPF). *Thromb Res.* 118(4): 463–9.
- [17] **Cannavo I et al. (2010):** Assessment of an immature platelet fraction (IPF%) in the diagnosis of thrombocytopenia. *Ann Biol Clin (Paris).* 68(4): 415–20.
- [18] **Strauss G et al. (2011):** Immature platelet count: a simple parameter for distinguishing thrombocytopenia in pediatric acute lymphocytic leukemia from immune thrombocytopenia. *Pediatr Blood Cancer.* 57(4): 641–7.
- [19] **Adly AA et al. (2015):** Evaluation of the immature platelet fraction in the diagnosis and prognosis of childhood immune thrombocytopenia. *Platelets.* 26(7): 645–50.
- [20] **Sakuragi M et al. (2015):** Clinical significance of IPF% or RP% measurement in distinguishing primary immune thrombocytopenia from aplastic thrombocytopenic disorders. *Int J Hematol.* 101(4): 369–75.
- [21] **Cremer M et al. (2016):** Thrombocytopenia and platelet transfusion in the neonate. *Semin Fetal Neonatal Med.* 21(1): 10–8.
- [22] **Miyazaki K et al. (2015):** Immature platelet fraction measurement is influenced by platelet size and is a useful parameter for discrimination of macrothrombocytopenia. *Hematology.* 20(10): 587–92.
- [23] **Sokolic R et al. (2015):** Assessment of immature platelet fraction in the diagnosis of Wiskott-Aldrich syndrome. *Front Pediatr.* 3: 1–10.
- [24] **Gerrits AJ et al. (2015):** Effects of eltrombopag on platelet count and platelet activation in Wiskott-Aldrich syndrome/X-linked thrombocytopenia. *Blood* 126(11): 1367–78.
- [25] **Tanaka Y et al. (2014):** Performance Evaluation of Platelet Counting by Novel Fluorescent Dye Staining in the XN-Series Automated Hematology Analyzers. *J Clin Lab Anal.* 28(5): 341–8.
- [26] **Park S et al. (2014):** The Sysmex XN-2000 Hematology Autoanalyzer Provides a Highly Accurate Platelet Count than the Former Sysmex XE-2100 System Based on Comparison with the CD41/CD61 Immunoplatelet Reference Method of Flow Cytometry. *Ann Lab Med.* 34(6): 471–4.

You can download our white papers from our website:  
[www.sysmex-europe.com/whitepapers](http://www.sysmex-europe.com/whitepapers)